


Product datasheet for **MO20021-100**

Cathepsin G (CTSG) Mouse Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Polyclonal
Formulation:	State: Supernatant
Gene Name:	cathepsin G
Database Link:	Entrez Gene 1511 Human
Synonyms:	CTSG, VSP

Product images:



Cathepsin G staining of Human Tonsil. Note intense membrane staining of polymorphonuclear leukocytes. Paraffin section.



Cathepsin G staining of Human Tonsil. Note intense membrane staining of polymorphonuclear leukocytes. Paraffin section.



[View online »](#)