

Product datasheet for **MO22104-100**

Neurofilament (NEFL) Mouse Monoclonal Antibody [Clone ID: DA2]

Product data:

Product Type:	Primary Antibodies
Clone Name:	DA2
Applications:	ELISA, IF, IHC, WB
Reactivity:	Chicken, Feline, Human, Mouse, Rat
Host:	Mouse
Isotype:	IgG
Clonality:	Monoclonal
Formulation:	State: Supernatant
Gene Name:	neurofilament, light polypeptide
Database Link:	Entrez Gene 4747 Human
Synonyms:	NF-L, NEFL, NF68, 68 kDa neurofilament protein



[View online »](#)