


Product datasheet for **MO22110-100**

Alpha Fodrin (SPTAN1) (676-1043) Mouse Monoclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IF, WB
Reactivity:	Human, Mouse, Rat
Host:	Mouse
Isotype:	IgG
Clonality:	Monoclonal
Formulation:	State: Aff - Purified
Gene Name:	Homo sapiens spectrin alpha, non-erythrocytic 1 (SPTAN1), transcript variant 1
Database Link:	Entrez Gene 20740 Mouse Entrez Gene 64159 Rat Entrez Gene 6709 Human
Synonyms:	Spectrin non-erythroid alpha chain, Alpha-II spectrin, Fodrin alpha chain, SPTA2
Protein Families:	Druggable Genome
Protein Pathways:	Tight junction

Product images:



[View online »](#)

