


## Product datasheet for **MO22941-100**

### Th Mouse Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	IF, IHC, WB
Reactivity:	Mouse, Primate, Rat
Host:	Mouse
Isotype:	IgG1
Clonality:	Polyclonal
Formulation:	State: Aff - Purified
Gene Name:	tyrosine hydroxylase
Database Link:	<a href="#">Entrez Gene 25085 Rat</a>
Synonyms:	Tyrosine 3-hydroxylase, TYH


#### Product images:



TH Staining in the Substantia Nigra.



[View online »](#)



TH Staining in the Substantia Nigra.