


Product datasheet for **MO25012-100**

14 3 3 gamma (YWHAG) Mouse Monoclonal Antibody [Clone ID: HS23]

Product data:

Product Type:	Primary Antibodies
Clone Name:	HS23
Applications:	WB
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Formulation:	State: Ascites
Gene Name:	Homo sapiens tyrosine 3-monooxygenase/tryptophan 5-monooxygenase activation protein gamma (YWHAG)
Database Link:	Entrez Gene 7532 Human
Synonyms:	YWHAG, KCIP-1
Protein Families:	Druggable Genome
Protein Pathways:	Cell cycle, Neurotrophin signaling pathway, Oocyte meiosis


Product images:



WB Analysis. HeLa whole cell lysate.



[View online »](#)



WB Analysis. HeLa whole cell lysate.