


## Product datasheet for MO25024-100

### Bestrophin (BEST1) Mouse Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	IF, WB
Reactivity:	Human, Porcine, Primate
Host:	Mouse
Isotype:	IgG
Clonality:	Polyclonal
Formulation:	State: Ascites
Gene Name:	bestrophin 1
Database Link:	<a href="#">Entrez Gene 7439 Human</a>
Synonyms:	BEST1, VMD2, TU15B

#### Product images:



[View online »](#)

