


## Product datasheet for MO25046-100

### EGLN1 (1-50) Mouse Monoclonal Antibody [Clone ID: 366G/76/3]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	366G/76/3
Applications:	IHC, WB
Reactivity:	Bovine, Human, Rat
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Formulation:	State: Aff - Purified
Gene Name:	Homo sapiens egl-9 family hypoxia-inducible factor 1 (EGLN1)
Database Link:	<a href="#">Entrez Gene 308913 Rat</a> <a href="#">Entrez Gene 54583 Human</a>
Synonyms:	EGLN-1, C1orf12, HIF-PH2, HIF-prolyl hydroxylase 2, HPH-2, SM-20
Protein Pathways:	Pathways in cancer, Renal cell carcinoma


#### Product images:



Western blot analysis in HeLa whole cell extracts using HIF Prolyl Hydroxylase 2 Antibody.



[View online »](#)



Western blot analysis in HeLa whole cell extracts using HIF Prolyl Hydroxylase 2 Antibody.