


Product datasheet for **MO50007-100**

Cacna1d (859-875) Mouse Monoclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IHC, IP, WB
Reactivity:	Human, Mouse, Rat
Host:	Mouse
Isotype:	IgG2a
Clonality:	Monoclonal
Formulation:	State: Supernatant
Gene Name:	calcium voltage-gated channel subunit alpha1 D
Database Link:	Entrez Gene 29716 Rat
Synonyms:	CACH3, CACN4, CACNL1A2, CCHL1A2, Cav1.3

Product images:



Cav1.3 Ca²⁺ channel staining of adult rat hippocampus.




[View online »](#)




CAV1.3 western blots of adult rat brain membrane.



Cav1.3Ca²⁺ channel western blots of adult rat brain membrane.



CAV1.3 staining of adult rat hippocampus.



Cav1.3 Ca²⁺ channel staining of adult rat hippocampus.



Cav1.3Ca²⁺ channel western blots of adult rat brain membrane.