

Product datasheet for **MO50015-100**

Kcna1 (458-476) Mouse Monoclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IF, IHC, IP, WB
Reactivity:	Human, Mouse, Rat
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Formulation:	State: Aff - Purified
Gene Name:	potassium voltage-gated channel subfamily A member 1
Database Link:	Entrez Gene 24520 Rat
Synonyms:	HUK1, HBK1

Product images:






Image: Kv1.1 staining of adult rat hippocampus.




[View online »](#)



Kv1.1 western blots of adult rat brain membranes (RBM), and adult mouse brain membranes from Kv1.1 knockout (MBM-Kv1.1-KO) and wild-type (MBM-WT) mice.



Kv1.1 western blots of adult rat brain membranes (RBM), and adult mouse brain membranes from Kv1.1 knockout (MBM-Kv1.1-KO) and wild-type (MBM-WT) mice.



Kv1.1 staining of adult rat hippocampus.





Image: Kv1.1 staining of adult rat hippocampus.



Kv1.1 western blots of adult rat brain membranes (RBM), and adult mouse brain membranes from Kv1.1 knockout (MBM-Kv1.1-KO) and wild-type (MBM-WT) mice.