

Product datasheet for **MP200058**

Htr4 Mouse qPCR Primer Pair (NM_008313)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	15562
Forward Sequence:	TGCTCACGTTCTTGCAGTGGT
Reverse Sequence:	GTCAGCAAAGGCGAGAGACACA
ACCN:	BC132189 , BC132191 , NM_008313 , NM_008313.1 , NM_008313.2 , NM_008313.3 , NM_008313.4 , BC145159 , BC145161 , BC148470 , BC153066
UniProt ID:	P97288
Synonyms:	5-HT-4; 5-HT4; 5-HT<4L>; 5HTR4
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)