

Product datasheet for **MP200071**

Htr1d Mouse qPCR Primer Pair (NM_008309)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	15552
Forward Sequence:	CTGGAGTACAGCAAGCGTCGAA
Reverse Sequence:	TGTTACCAGGCAGTCGGACAT
ACCN:	BC103532 , BC103535 , BC103536 , NM_008309 , NM_008309.1 , NM_008309.2 , NM_008309.3 , NM_008309.4 , NM_008309.5 , BB832969 , BY728591
UniProt ID:	Q61224
Synonyms:	5-HT-1D; 5-HT1D; AI853647; Gpcr14; Htr1db
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)