

Product datasheet for **MP200088**

Nt5c2 Mouse qPCR Primer Pair (NM_029810)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	76952
Forward Sequence:	GTGGCAGTCCTACTTTGACCTG
Reverse Sequence:	CCTGAGTAGACAATGCCGTGCT
ACCN:	<u>BC064760</u> , <u>NM_029810</u> , <u>NM_029810.1</u> , <u>NM_029810.2</u> , <u>NM_029810.3</u> , <u>NM_029810.4</u>
UniProt ID:	<u>Q3V1L4</u>
Synonyms:	2010002123Rik; cN-II; CnII; ENSMUSG00000025039; Gm9751; Gmp; Nt5b; Pnt5
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)