

Product datasheet for **MP200089**

Nt5c3 Mouse qPCR Primer Pair (NM_026004)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	107569
Forward Sequence:	GTTGATCCTGTTCTCACCGTGG
Reverse Sequence:	CGTCAGAGTCAGCCACAATCTC
ACCN:	<u>BC038029</u> , <u>NM_026004</u> , <u>NM_026004.1</u> , <u>NM_026004.2</u> , <u>NM_026004.3</u> , <u>BX527194</u> , <u>BY080680</u>
UniProt ID:	<u>Q9D020</u>
Synonyms:	1600024P05Rik; 2610206B05Rik; 3110004A18Rik; AI842776; cN-III; lupin; Nt5c3a; Nt5y; p36
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)