

Product datasheet for **MP200159**

Akap8 Mouse qPCR Primer Pair (NM_019774)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	56399
Forward Sequence:	TGAGAAACGGCGTCAGGAGTTG
Reverse Sequence:	CGCAGGAATCAGCATGTCACAG
ACCN:	<u>BC125521</u> , <u>NM_019774</u> , <u>NM_019774.1</u> , <u>NM_019774.2</u> , <u>NM_019774.3</u> , <u>NM_019774.4</u> , <u>NM_019774.5</u> , <u>BC042602</u> , <u>BC054433</u> , <u>BC060170</u> , <u>BC060209</u> , <u>BC078448</u> , <u>BC137746</u>
UniProt ID:	<u>Q9DBR0</u>
Synonyms:	1200016A02Rik; AA673585; AKA; AKAP-8; AKAP95; AU015639
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)