

Product datasheet for **MP200177**

Abhd6 Mouse qPCR Primer Pair (NM_025341)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	66082
Forward Sequence:	GGCATTGTGTCGCTTTCTCTCC
Reverse Sequence:	CGGAAGGAGTAACAGAACTGATAG
ACCN:	BC027011 , NM_025341 , NM_025341.1 , NM_025341.2 , NM_025341.3 , NM_025341.4 , BC051086 , BE624233 , BQ176347 , BY253400
UniProt ID:	Q8R2Y0
Synonyms:	0610041D24Rik; AA673485; AV065425
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)