

Product datasheet for **MP200189**

Acat1 Mouse qPCR Primer Pair (NM_144784)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	110446
Forward Sequence:	GCAGGGAAGTTTGCCAGTGAGA
Reverse Sequence:	GAACACGGTCTTGAGCTTTGGC
ACCN:	<u>BC024763</u> , <u>NM_144784</u> , <u>NM_144784.1</u> , <u>NM_144784.2</u> , <u>NM_144784.3</u> , <u>BC024429</u> , <u>BC029345</u>
UniProt ID:	<u>Q8QZT1</u>
Synonyms:	6330585C21Rik; Acat
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)