

## Product datasheet for **MP200203**

### Acp5 Mouse qPCR Primer Pair (NM\_007388)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	11433
Forward Sequence:	GCGACCATTGTTAGCCACATACG
Reverse Sequence:	CGTTGATGTCGCACAGAGGGAT
ACCN:	<a href="#">BC012911</a> , <a href="#">BC019160</a> , <a href="#">BC029644</a> , <a href="#">NM_007388</a> , <a href="#">NM_007388.1</a> , <a href="#">NM_007388.2</a> , <a href="#">NM_007388.3</a> , <a href="#">BG076865</a> , <a href="#">BQ031772</a> , <a href="#">BY084869</a> , <a href="#">BY107232</a>
UniProt ID:	<a href="#">Q05117</a>
Synonyms:	TRACP; TRAP
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)