

Product datasheet for **MP200275**

Acsf3 Mouse qPCR Primer Pair (NM_144932)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	257633
Forward Sequence:	CACTTCCGATGGCTGGTTCAGA
Reverse Sequence:	GCGCTGACTTTATAGCCTCCAG
ACCN:	<u>BC022709</u> , <u>NM_144932</u> , <u>NM_144932.1</u> , <u>NM_144932.2</u> , <u>NM_144932.3</u> , <u>NM_144932.4</u> , <u>BC021611</u> , <u>BC156691</u> , <u>BG920359</u> , <u>BI135241</u> , <u>BQ177632</u>
UniProt ID:	<u>Q3URE1</u>
Synonyms:	BB101783
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)