

Product datasheet for **MP200298**

Acot7 Mouse qPCR Primer Pair (NM_133348)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	70025
Forward Sequence:	GGACACCTGTTCTCAGATTCCG
Reverse Sequence:	CCAGCCACATTGGCATCATCTG
ACCN:	<u>BC013507</u> , <u>NM_133348</u> , <u>NM_133348.1</u> , <u>NM_133348.2</u> , <u>BG069979</u> , <u>BQ748599</u> , <u>BY134190</u> , <u>BY327697</u>
UniProt ID:	<u>Q91V12</u>
Synonyms:	2410041A17Rik; Ach1; Act; Bach; Cte-II; CTE-IIa; Lach1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)