

Product datasheet for **MP200308**

Acad9 Mouse qPCR Primer Pair (NM_172678)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	229211
Forward Sequence:	CTCCAACACGATGTATGCTCGG
Reverse Sequence:	GCCTTCTGTTCTCATTGCCGA
ACCN:	<u>BC031137</u> , <u>BC032213</u> , <u>BC033277</u> , <u>NM_172678</u> , <u>NM_172678.1</u> , <u>NM_172678.2</u> , <u>NM_172678.3</u> , <u>NM_172678.4</u> , <u>BC004737</u>
UniProt ID:	<u>Q8JZN5</u>
Synonyms:	2600017P15Rik; 4732402K02; C630012L17Rik; NPD002
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)