

Product datasheet for **MP200331**

Ap1b1 Mouse qPCR Primer Pair (NM_007454)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	11764
Forward Sequence:	GGTTCAGGAAGCCATTGTGGTC
Reverse Sequence:	TCAGGCTCATCCAGAGAGTCCA
ACCN:	BC008513 , NM_007454 , NM_007454.1 , NM_007454.2 , BC066827 , BP765289 , BY253407
UniProt ID:	O35643
Synonyms:	Adtb1; b2b1660Clo
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)