

Product datasheet for **MP200354**

Ap4e1 Mouse qPCR Primer Pair (NM_175550)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	108011
Forward Sequence:	TCTGGAAACAGTGCCGAGACAG
Reverse Sequence:	CCAGTTCCACTCTCCTTCTTGG
ACCN:	<u>NM_175550</u> , <u>NM_175550.1</u> , <u>NM_175550.2</u> , <u>NM_175550.3</u> , <u>BC042530</u> , <u>BC152829</u> , <u>BY037321</u>
UniProt ID:	<u>Q80V94</u>
Synonyms:	2310033A20Rik; 9930028M04Rik; AV087807
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

