

Product datasheet for **MP200397**

Ak4 Mouse qPCR Primer Pair (NM_009647)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	11639
Forward Sequence:	GAAGCAGTTGCTGCCAGGCTAA
Reverse Sequence:	GCCAGATTCTGTTAGTCTCCGTC
ACCN:	<u>BC086663</u> , <u>NM_009647</u> , <u>NM_009647.1</u> , <u>NM_009647.2</u> , <u>NM_009647.3</u> , <u>NM_009647.4</u> , <u>NM_009647.5</u> , <u>BC086791</u> , <u>BY203078</u> , <u>BY718404</u> , <u>BY740484</u>
UniProt ID:	<u>Q9WUR9</u>
Synonyms:	Ak-3; Ak-4; Ak3; Ak3l1; AK 4; D4Ertd274e
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)