

Product datasheet for **MP200502**

Aarsd1 Mouse qPCR Primer Pair (NM_144829)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	69684
Forward Sequence:	TGTGGTTCTGGAAGACACGCTG
Reverse Sequence:	GACAGAGGTGACTCCGTGAAGT
ACCN:	<u>BC005711</u> , <u>NM_144829</u> , <u>NM_144829.1</u>
UniProt ID:	<u>Q3THG9</u>
Synonyms:	1110069E20Rik; 2310044P18Rik; AA589600; AlaX; AlaXp-II
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 μ l dH ₂ O to make a final concentration of 10 μ M.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)