

## Product datasheet for **MP200511**

### Aldh1l1 Mouse qPCR Primer Pair (NM\_027406)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	107747
Forward Sequence:	CTTCATAGGCGGCGAGTTTGTG
Reverse Sequence:	CGCCTTGTC AACATCACTCACC
ACCN:	<u><a href="#">BC024055</a></u> , <u><a href="#">BC025939</a></u> , <u><a href="#">BC028817</a></u> , <u><a href="#">BC030722</a></u> , <u><a href="#">BC030723</a></u> , <u><a href="#">BC030727</a></u> , <u><a href="#">BC030730</a></u> , <u><a href="#">NM_027406</a></u> , <u><a href="#">NM_027406.1</a></u> , <u><a href="#">BC030722.1</a></u> , <u><a href="#">BC022112</a></u> , <u><a href="#">NM_027406.2</a></u>
UniProt ID:	<u><a href="#">Q8R0Y6</a></u>
Synonyms:	1810048F20Rik; FDH; Fthfd; Neut2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)