

Product datasheet for **MP200515**

Aldh2 Mouse qPCR Primer Pair (NM_009656)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	11669
Forward Sequence:	GCTGTTGTACCGATTGGCGGAT
Reverse Sequence:	GCGGAGACATTTTCAGGACCATG
ACCN:	<u>BC005476, NM_009656, NM_009656.1, NM_009656.2, NM_009656.3, NM_009656.4</u>
UniProt ID:	<u>P47738</u>
Synonyms:	Ahd-5; AHD-M1; Ahd5; ALDH-E2; ALDHI
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)