

Product datasheet for **MP200532**

Aldh5a1 Mouse qPCR Primer Pair (NM_172532)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	214579
Forward Sequence:	CCTTTACTGGCTCAACAGCAACG
Reverse Sequence:	CTGGTCTACGTTGGCACTGTCA
ACCN:	<u>NM_172532</u> , <u>NM_172532.1</u> , <u>NM_172532.2</u> , <u>NM_172532.3</u> , <u>BC138703</u> , <u>BC030390</u>
UniProt ID:	<u>Q8BWF0</u>
Synonyms:	6330403E24Rik; Ahd-1; Ahd1; D630032B01Rik; OTTMUSG00000000613; SSADH; Ssdh1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)