

Product datasheet for **MP200589**

Asic2 Mouse qPCR Primer Pair (NM_007384)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	11418
Forward Sequence:	GGAGATCATGCTGGACATTCAGC
Reverse Sequence:	GCTCTTGGATGAAAGGTGGCTC
ACCN:	<u>NM_007384</u> , <u>NM_007384.1</u> , <u>NM_007384.2</u> , <u>NM_007384.3</u> , <u>BC138529</u> , <u>BB342826</u> , <u>BC038551</u> , <u>BC138530</u>
UniProt ID:	<u>Q925H0</u>
Synonyms:	Accn1; ACIC2; BNaC1; BNaC1a; BNC1; Mdeg
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)