

Product datasheet for **MP200701**

Anks1 Mouse qPCR Primer Pair (NM_181413)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	224650
Forward Sequence:	TCCTCTCCAACACAGACGCCAT
Reverse Sequence:	GAATGGCAGGACACAGCGTTGT
ACCN:	BC050847 , NM_181413 , NM_181413.1 , NM_181413.2 , NM_181413.3 , NM_181413.4 , BC010633 , BC027740 , BC057331
UniProt ID:	P59672
Synonyms:	Anks1a; mKIAA0229
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)