

Product datasheet for **MP200714**

Asb4 Mouse qPCR Primer Pair (NM_023048)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	65255
Forward Sequence:	TCCATCCTCTGTGCCAAGCAGT
Reverse Sequence:	CGTAGAAGGCAACCAGCTCAGA
ACCN:	<u>BC046819</u> , <u>BC054745</u> , <u>BC056440</u> , <u>NM_023048</u> , <u>NM_023048.1</u> , <u>NM_023048.2</u> , <u>NM_023048.3</u> , <u>NM_023048.4</u> , <u>NM_023048.5</u>
UniProt ID:	<u>Q9WV71</u>
Synonyms:	8430401O13Rik; AV117194; AW060255
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)