

Product datasheet for **MP200813**

Apoc3 Mouse qPCR Primer Pair (NM_023114)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	11814
Forward Sequence:	AAGACGGTCCAGGATGCGCTAA
Reverse Sequence:	GTTGGTCTCAGGTTAGAATCC
ACCN:	<u>BC021776, NM_023114, NM_023114.1, NM_023114.2, NM_023114.3, NM_023114.4</u>
UniProt ID:	<u>P33622</u>
Synonyms:	apo-CIII; apoC-III
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)