

Product datasheet for **MP200815**

Apod Mouse qPCR Primer Pair (NM_007470)

Product data:

| | |
|-------------------|---|
| Product Type: | qPCR Primer Pairs |
| Gene ID: | 11815 |
| Forward Sequence: | GGTGAAGCCAAACAGAGCAACG |
| Reverse Sequence: | CAGGAGTACACGAGGGCATAGT |
| ACCN: | <u>NM_007470</u> , <u>NM_007470.1</u> , <u>NM_007470.2</u> , <u>NM_007470.3</u> , <u>NM_007470.4</u> , <u>BC145907</u> , <u>BC145909</u> , <u>BY596441</u> |
| UniProt ID: | <u>P51910</u> |
| Components: | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM. |
| Quality Control: | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage: | Store at -20°C. |
| Stability: | The primer mix is stable for one year from date of shipping. |



[View online »](#)