

Product datasheet for **MP200946**

Arsk Mouse qPCR Primer Pair (NM_029847)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	77041
Forward Sequence:	GTGTGCTGAGACAGATGCCATG
Reverse Sequence:	AAACTGGCGGTGTTCCATAGCC
ACCN:	BC046790 , BC058351 , NM_029847 , NM_029847.1 , NM_029847.2 , NM_029847.3 , NM_029847.4 , BN000751
UniProt ID:	Q9D2L1
Synonyms:	2810429K17Rik; 4833414G15Rik
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)