

Product datasheet for **MP200961**

Alg13 Mouse qPCR Primer Pair (NM_026247)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	67574
Forward Sequence:	TGGTGCCCAAACCATTCGGTAC
Reverse Sequence:	CAGACTCTCCAAACAGCTTCCTG
ACCN:	<u>NM_026247</u> , <u>NM_026247.1</u> , <u>NM_026247.2</u> , <u>NM_026247.3</u> , <u>BC137692</u> , <u>BC053004</u> , <u>BC059024</u>
UniProt ID:	<u>Q9D8C3</u>
Synonyms:	2810046O15Rik; 4833435D08Rik; Glt28d1; MDS031; OTTMUSG00000018917
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)