

Product datasheet for **MP200978**

Aga Mouse qPCR Primer Pair (NM_001005847)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	11593
Forward Sequence:	GAAACCACCTGGATGCCATGA
Reverse Sequence:	GTGTGGTATGCTCCAGGACTCT
ACCN:	BC132227 , BC132229 , NM_001005847 , NM_001005847.1 , NM_001005847.2 , BC048230 , BC145297 , BE656262 , BQ032917 , BQ571087
UniProt ID:	Q64191
Synonyms:	AW060726
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)