

## Product datasheet for **MP200979**

### Dars Mouse qPCR Primer Pair (NM\_177445)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	226414
Forward Sequence:	GCATCTGCCATCTCTCCGAGA
Reverse Sequence:	TGTGGAGACTGAGCAAGGTAGG
ACCN:	<u><a href="#">BC008638</a></u> , <u><a href="#">NM_177445</a></u> , <u><a href="#">NM_177445.1</a></u> , <u><a href="#">NM_177445.2</a></u> , <u><a href="#">NM_177445.3</a></u> , <u><a href="#">NM_177445.4</a></u> , <u><a href="#">NM_177445.5</a></u> , <u><a href="#">BB576692</a></u> , <u><a href="#">BY276340</a></u> , <u><a href="#">BY308833</a></u>
UniProt ID:	<u><a href="#">Q922B2</a></u>
Synonyms:	5730439G15Rik
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)