

Product datasheet for **MP200999**

Atxn1 Mouse qPCR Primer Pair (NM_009124)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	20238
Forward Sequence:	ACAGTTTGCTGTTGGTGAACACC
Reverse Sequence:	ACACGGCAGATCAAAGAGCTGG
ACCN:	<u>BC058178</u> , <u>NM_009124</u> , <u>NM_009124.1</u> , <u>NM_009124.2</u> , <u>NM_009124.3</u> , <u>NM_009124.5</u> , <u>NM_009124.6</u> , <u>BE993799</u>
UniProt ID:	<u>P54254</u>
Synonyms:	2900016G23Rik; Atx1; C85907; ENSMUSG00000074917; Gm10786; Scal
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

