

Product datasheet for **MP201023**

Acly Mouse qPCR Primer Pair (NM_134037)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	104112
Forward Sequence:	AGGAAGTGCCACCTCCAACAGT
Reverse Sequence:	CGCTCATCACAGATGCTGGTCA
ACCN:	BC056378 , NM_134037 , NM_134037.1 , NM_134037.2 , NM_134037.3 , BC005533 , BC021502 , BC065805 , BE623317
UniProt ID:	Q91V92
Synonyms:	A730098H14Rik; AW538652
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)