

Product datasheet for **MP201037**

Atp5d Mouse qPCR Primer Pair (NM_025313)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	66043
Forward Sequence:	ACTGGAGCCTTTGGCATCTTG
Reverse Sequence:	AGTCGGCATTACAGTGACGGA
ACCN:	<u>BC008273</u> , <u>NM_025313</u> , <u>NM_025313.1</u> , <u>NM_025313.2</u> , <u>BE951291</u> , <u>BM948933</u> , <u>BY451544</u>
UniProt ID:	<u>Q9D3D9</u>
Synonyms:	0610008F14Rik; 1500000I11Rik; AA960090; AI876556; AU020773; C85518
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)