

Product datasheet for **MP201050**

Atad3a Mouse qPCR Primer Pair (NM_179203)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	108888
Forward Sequence:	GTGACTGCCATGCACAAGGTCT
Reverse Sequence:	AGTCGCTCGTTTCCTGAGGAAG
ACCN:	<u>BC023301</u> , <u>BC058373</u> , <u>BC060036</u> , <u>NM_179203</u> , <u>NM_179203.1</u> , <u>NM_179203.2</u> , <u>NM_179203.3</u> , <u>BB839884</u> , <u>BC047047</u> , <u>BF318525</u>
UniProt ID:	<u>Q92511</u>
Synonyms:	2400004H09Rik; Atad3; mKIAA1273; Tob3
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)