

Product datasheet for **MP201054**

Atp13a1 Mouse qPCR Primer Pair (NM_133224)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	170759
Forward Sequence:	CCTCAAAGCTGACTCCAAGGCT
Reverse Sequence:	GCAAGATGAGCGTGTGGGCTTT
ACCN:	<u>NM_133224</u> , <u>NM_133224.1</u> , <u>NM_133224.2</u> , <u>BC138721</u> , <u>BC022689</u> , <u>BC049817</u> , <u>BC078635</u> , <u>BC138722</u>
UniProt ID:	<u>Q9EPE9</u>
Synonyms:	Atp13a; catp; Cgi152
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.


[View online »](#)