

Product datasheet for **MP201103**

Atp6v1g2 Mouse qPCR Primer Pair (NM_023179)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	66237
Forward Sequence:	GGAACAGGAGTTTCAGAGCAAGC
Reverse Sequence:	TTCCTCTGCTGGAACTCTGCA
ACCN:	<u>BC020190</u> , <u>BC062380</u> , <u>NM_023179</u> , <u>NM_023179.1</u> , <u>NM_023179.2</u> , <u>NM_023179.3</u> , <u>BC025190</u> , <u>BC049587</u> , <u>BU558667</u>
UniProt ID:	<u>Q9WTT4</u>
Synonyms:	1500002D01Rik; Atp6g2; NG38; VAG2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)