

Product datasheet for **MP201123**

Abca4 Mouse qPCR Primer Pair (NM_007378)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	11304
Forward Sequence:	TATGCCAGCCTTTTCCGAGAGC
Reverse Sequence:	TCCAGCATCCTCTGTGACCTTC
ACCN:	BC057853 , NM_007378 , NM_007378.1 , BC043937
UniProt ID:	O35600
Synonyms:	Abc1; Abc10; Abcr; AW050280; D430003I15Rik; R; RmP
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)