

Product datasheet for **MP201133**

Abcb1b Mouse qPCR Primer Pair (NM_011075)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	18669
Forward Sequence:	GGTGGTGTCATTGTGGAGCAAG
Reverse Sequence:	GCATCAGTGTCACCTCTGGGATC
ACCN:	<u>NM_011075</u> , <u>NM_011075.1</u> , <u>NM_011075.2</u> , <u>BC141363</u> , <u>BC150811</u> , <u>BG080311</u> , <u>BQ829721</u> , <u>BY435708</u>
UniProt ID:	<u>P06795</u>
Synonyms:	Abcb1; mdr; Mdr1; Mdr1b; Pgy-1; Pgy1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)