

Product datasheet for **MP201142**

Abcc3 Mouse qPCR Primer Pair (NM_029600)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	76408
Forward Sequence:	ACTTCCTCCGAAACTACGCACC
Reverse Sequence:	GCTGGCTCATTGTCTGTCAAGT
ACCN:	<u>NM_029600</u> , <u>NM_029600.1</u> , <u>NM_029600.2</u> , <u>NM_029600.3</u> , <u>BC150788</u> , <u>BC043728</u> , <u>BC046560</u> , <u>BC048825</u> , <u>BC058185</u> , <u>NM_029600.4</u>
UniProt ID:	<u>B2RX12</u>
Synonyms:	1700019L09Rik; ABC31; MLP2; MOAT-D; MRP3
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)