

Product datasheet for **MP201157**

Abcg2 Mouse qPCR Primer Pair (NM_011920)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	26357
Forward Sequence:	CAGTTCTCAGCAGCTCTTCGAC
Reverse Sequence:	TCCTCCAGAGATGCCACGGATA
ACCN:	<u>BC053730, NM_011920, NM_011920.1, NM_011920.2, NM_011920.3</u>
UniProt ID:	<u>Q7TMS5</u>
Synonyms:	ABC15; ABCP; AI428558; BCRP; Bcrp1; MXR; MXR1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)