

Product datasheet for **MP201164**

Atrnl1 Mouse qPCR Primer Pair (NM_181415)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	226255
Forward Sequence:	GCAGGAATGTCAGGCTTTGCCA
Reverse Sequence:	GAGGCGAGATTTCCACTGAGGT
ACCN:	BC050020 , NM_181415 , NM_181415.1 , NM_181415.2 , NM_181415.3 , NM_181415.4 , BC022118 , BC027764 , BC030872 , BM899055 , BM944131 , BU610577
UniProt ID:	Q6A051
Synonyms:	AI504415; Alp; Atrnl; Atrnlp; AW555641
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.


[View online »](#)