

Product datasheet for **MP201169**

Aurkb Mouse qPCR Primer Pair (NM_011496)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	20877
Forward Sequence:	CTTCTACGACCAGCAGAGGATC
Reverse Sequence:	GGCATCTGACAGTTCCTCCATG
ACCN:	<u>BC003261</u> , <u>NM_011496</u> , <u>NM_011496.1</u> , <u>NM_011496.2</u> , <u>BY005259</u>
UniProt ID:	<u>O70126</u>
Synonyms:	Al; Aik2; AIM-1; Aim1; AIRK2; AL022959; Ark2; AurB; I; IPL1; Stk; STK-; STK-1; Stk1; Stk5; Stk12
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)