

Product datasheet for **MP201185**

Atg9a Mouse qPCR Primer Pair (NM_001003917)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	245860
Forward Sequence:	GTTAGCTGTGGAACACGTCCTC
Reverse Sequence:	GCAAGAATCACTCGGAGCAGCT
ACCN:	<u>BC079884</u> , <u>NM_001003917</u> , <u>NM_001003917.1</u> , <u>NM_001003917.2</u> , <u>NM_001003917.3</u> , <u>NM_001003917.4</u> , <u>BC034139</u> , <u>BC056488</u> , <u>BC082307</u> , <u>BK004020</u>
UniProt ID:	<u>Q68FE2</u>
Synonyms:	Apg9l1; Atg9; Atg9l1; AU019532
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60 °C for 1 min; Stage 4: Melting curve: 95 °C for 15 sec, 60 °C for 15 sec, 95 °C for 15 sec.
Storage:	Store at -20 °C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)