

Product datasheet for **MP201212**

Bbs4 Mouse qPCR Primer Pair (NM_175325)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	102774
Forward Sequence:	GCCTGTGCTATTCCAGAAAGTCC
Reverse Sequence:	TCTTCCAGTCGAAGGGTGCCAA
ACCN:	<u>NM_175325</u> , <u>NM_175325.1</u> , <u>NM_175325.2</u> , <u>NM_175325.3</u> , <u>BC145771</u> , <u>BC055797</u> , <u>BC089507</u> , <u>BC092531</u> , <u>BC145773</u>
UniProt ID:	<u>Q8C1Z7</u>
Synonyms:	AW537059; AW742241; D9Ert464e
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)